

SAM

Introducing the next level in patient care, rest and rehabilitation. The Sam Patient Chair is available in both standard and headrest models with either a flex action or a full motion recline system.



Model	Description	Dimensions				
		Overall			Seat	
		H	W	D	W	D
101M	Patient Chair Supporting Flex Action Sitting Comfort. Replaceable Upholstered Seat and Back. Pivoting Soft Urethane Armrests. Mobilizes On Casters When Tilted.	41	30	29	23	19
102M	Patient Chair with Headrest Supporting Flex Action Sitting Comfort. Replaceable Upholstered Seat and Back. Pivoting Soft Urethane Armrests. Mobilizes On Casters When Tilted.	48	30	30	23	19
103M	Patient Chair Full Motion Recline System w/ Pneumatic Seat Lift. Replaceable Full Upholstered Seat and Back. Pivoting Soft Urethane Armrests. Mobilizes on Casters When Tilted.	41	30	29	23	19
104M	Patient Chair with Headrest Full Motion Recline System w/ Pneumatic Seat Lift. Replaceable Full Upholstered Seat and Back. Pivoting Soft Urethane Armrests. Mobilizes on Casters When Tilted.	48	30	30	23	19



Frame Finish: Silver or Mahogany
Arm Rests: Silver Grey or Mahogany

Features: * Flex-Action Suspension System (models 101-2)
 * Full Motion Recline System (models 103-4)
 * Pivoting Soft Urethane Arm Rests
 * Pneumatic Seat Lift (models 103-4)
 * Light and Compact For Handling Ease
 * Mobilizes On Casters When Tilted
 * Field Replaceable/Cleanable Components

List Prices - Cover Grades

COM

H	I	K	M	P	S	V	Y	L-1
1295	1314	1352	1390	1447	1504	1561	1618	1695

Wgt. **Yds.** ***Cubes** **Load**
 47 1.5 18 350

**Note: Cubes based on back shipped KD (4 bolts to install)*

1395	1421	1473	1525	1603	1681	1759	1837	1945
------	------	------	------	------	------	------	------	------

Wgt. **Yds.** ***Cubes** **Load**
 50 2 18 350

**Note: Cubes based on back shipped KD (4 bolts to install)*

1680	1699	1737	1775	1832	1889	1946	2003	2080
------	------	------	------	------	------	------	------	------

Wgt. **Yds.** ***Cubes** **Load**
 54 1.5 18 500

**Note: Cubes based on back shipped KD (4 bolts to install)*

1780	1806	1855	1910	1988	2066	2144	2222	2330
------	------	------	------	------	------	------	------	------

Wgt. **Yds.** ***Cubes** **Load**
 57 2 18 500

**Note: Cubes based on back shipped KD (4 bolts to install)*

